

TOE/TEE kit (Adult)

REF 1251-02

Description

TOE/TEE kit sterile.
A transducer cover (round tip) on an applicator, an extended 10 ml syringe filled with ultrasound transmission gel, a twist-lock and an adult bite guard.

Packaging information

| | |
|--------------|---------------|
| Packaging | Peel pouch |
| Selling unit | 02x24 kits |
| EAN13 (x24) | 8717371750589 |
| EAN13 (x48) | 8717371750596 |

Specifications

| | |
|--------------|--------------|
| Sterile | Yes (gamma) |
| Material | Polyurethane |
| Latex free | Yes |
| Tip Ø | 15 mm |
| Tapered to Ø | 32 mm |
| Total length | 1100 mm |



Probetection™ TOE/TEE kit.

Specially designed for Adult Trans Oesophageal /Trans Esophageal Echocardiography. (TOE/TEE)

The Probetection™ TOE/TEE is available as a complete, ready to use kit.

Palmedic offers a complete line of ultrasound transducer covers to protect the patients and healthcare workers against cross contamination.

Probetection™ features:

- Optimal prevention against cross contamination for patients and medical staff.
- Ready to use, including ultrasound gel and bite guard.
- Time saving for medical staff.
- Reduces probe down-time at HLD.
- Probe lasts longer because HLD is not needed between each patient.

Probetection™

Specially designed for Adult Trans Oesophageal or Trans Esophageal Echocardiography (TOE/TEE).

The Probetection™ TOE/TEE is available as a complete, ready to use kit.

Palmedic offers a complete line of ultrasound transducer covers to protect the patients and healthcare workers against cross contamination.

All TOE/TEE covers are made out of a latex free polyurethane material. All our TOE/TEE covers will fit on every TEE transducer available in the market.



Probetection™

TOE/TEE transducer kits and covers



Probetection™ kit

For which transducer?

The 1250/51 range is suitable for all ADULT TOE/TEE transducer brands used.

Amongst others the:

- Philips X7-2t live 3D TEE
- Philips S7-2 Omni TEE
- Esaote TEE022 Multiplane
- Esaote ST2612
- Esaote TEE122 Multiplane
- GE 4D TEE
- GE 6T TEE
- GE 6T-RS
- GE 6TC
- Siemens V5Ms

- Hitachi Aloka UST-5293-5
- Hitachi Aloka UST-52116
- Hitachi Aloka UST-52104
- Samsung Medison MPT4-7
- Toshiba TEE probe PEF-510MA
- Toshiba TEE probe PEF-510MB
- Toshiba TEE probe PEF-512MC
- Shison TEE
- SonoScape MPTEE
- SonoSite TEE
- SonoSite TEEEx
- Samsung Medison MPT4-7
- ZONARE P8-3TEE



Probetection™ bite guard